



PRO  
**BODYTITE**  
*by* INMODE



MINIMALLY INVASIVE RADIO-FREQUENCY  
ASSISTED LIPOLYSIS (RFAL)

---

*Surgical Results Without the Large Scars*

# BODYtite + FACEtite

FIVE TIMES FASTER THAN ANY OTHER THERMAL SYSTEM

BodyTite provides superior internal and external lipolysis, delivering body contouring results previously only achieved through more extensive surgical procedures. Powered by RFAL (Radio-Frequency Assisted Lipolysis), BodyTite provides unprecedented three-dimensional contraction of the dermis, subdermal connective tissue matrix and deep fat which is far superior than any other technology on the market.


BodyTite works by delivering rapid and uniform heating of the entire soft tissue matrix. Built-in safeguards including real time measurements of skin temperature, impedance monitoring, power cut-off and audible feedback help ensure patient safety.

Internal Temperature: | Operator adjustable between 50°C - 70°C

External Temperature: | Operator adjustable up to 42°C






	FACEtite	BODYtite	BODYTITE <sup>PRO</sup>
Power:	22 W	22 W	40 W
Cannula Length:	10 cm	17 cm	17 cm
Cannula Diameter:	1.8 mm	3.9 mm	3.9 mm
Depth of Treatment:	up to 25 mm (1")	up to 50 mm (2")	up to 50 mm (2")

RFAL Cannulas



# BODYtite + FACEtite

IS SAFE & DELIVERS TARGETED RESULTS

Functionality	Benefits	
Increased power to 40W	BodyTitePRO heats 2X faster than BodyTite and 5X faster than any other thermal system for total soft tissue contraction; not just for the skin, but the entire fibroseptal network	
Uniform thermal field provides an even distribution of heat	Controlled heating at critical temperatures allows contraction at 36.4-47%, far higher than any other thermal energy system <sup>1</sup>	
Immediate real-time audio & visual temperature feedback	Temperature feedback measures at 1000X per second to provide real time results, even in deep fat areas where laser temperature readers cannot reach	
100% energy concentration in the treatment zone	Allows for focused treatment of areas where operators want to avoid heating peripheral nerves or structures	
Physician pre-sets and controls heat gradient enabling optimal treatment results	Superficial contouring up to 25 mm (1 inch) and deep fat 50+mm (2 inches) allows treatment on various body areas	



## BENEFITS TO YOUR PRACTICE

- » BodyTite now offers twice as much power to improve treatment times
- » Engage patients with an innovative suite of minimally invasive procedures that deliver visible results with limited downtime
- » FaceTite delivers dramatic improvement in skin laxity (without the scars), which was previously only attainable with invasive procedures such as a facelift

## A.C.E. TECHNOLOGY

A contained thermal field ensures a controlled and safe treatment through A.C.E. (Acquire, Control, Extend) technology.

**ACQUIRE:** operators can acquire critical thermal temperatures that are best suited for contraction.

**CONTROL:** the operator programmable internal and external limit allows the treatment to be cut-off when the temperature is achieved.

**EXTEND:** maximize results by exposing the skin and adipose tissue to clinically therapeutic temperatures for long periods without overheating.

1. Chia, C. T., Neinstein, R. M., & Theodorou, S. J. (2017). Evidence-Based Medicine: Liposuction. *Plastic and Reconstructive Surgery*, 139(1). [http://journals.lww.com/plasreconsurg/Fulltext/2017/01000/Evidence\\_Based\\_Medicine\\_\\_\\_Liposuction.51.aspx](http://journals.lww.com/plasreconsurg/Fulltext/2017/01000/Evidence_Based_Medicine___Liposuction.51.aspx)

Duncan, D. I. (2013). Nonexcisional Tissue Tightening Creating Skin Surface Area Reduction During Abdominal Liposuction by Adding Radiofrequency Heating. *Aesthetic Surgery Journal*, 33(8), 1154-1166. <https://doi.org/10.1177/1090820X13505862>



# BODYTITE & FACETITE; LEADERS IN CLINICAL EXCELLENCE

**LINEAR CONTRACTION OBSERVED AT 6-MONTH FOLLOW-UP WAS MUCH MORE SIGNIFICANT THAN REPORTED WITH ANY OTHER TECHNOLOGY AND VARIED FROM 12.7% UP TO 47%, DEPENDING ON PATIENT AND TREATMENT VARIABLES.**

Christopher T. Chia, M.D., Ryan M. Neinstein, M.D. and Spero J. Theodorou, M.D.  
Plastic and Reconstructive Surgery, Evidence-Based Medicine: Liposuction

**THE MEAN SKIN SURFACE AREA REDUCTION WAS 25.8% IN REGIONS TREATED WITH RADIO-FREQUENCY PLUS SAL AT 6 WEEKS AND 36.4% AT 1 YEAR.**

D. Duncan, M.D., Aesthetic Surgery Journal, Nonexcisional Tissue Tightening: Creating Skin Surface Area Reduction During Abdominal Liposuction by Adding Radio-frequency Heating

For additional BodyTite and FaceTite clinical data visit: [www.inmodemd.com](http://www.inmodemd.com)

## JOIN THE REVOLUTION

Just some of the many plastic surgeons successfully integrating BodyTite into their patient solutions.

- |                         |                        |                           |                        |
|-------------------------|------------------------|---------------------------|------------------------|
| Min S. Ahn, MD          | John Q. Cook, MD       | Dennis Hurwitz, MD        | Nicholas Nikolov, MD   |
| Jason Altman, MD        | Frederic Corbin, MD    | Christopher Khorsandi, MD | Jason Pozner, MD       |
| Mark Anton, MD          | Mark Craig, MD         | Kenneth Kim, MD           | Leif Rogers, MD        |
| Robert Applebaum, MD    | Daniel Del Vecchio, MD | George Kouris, MD         | Pier Paolo Rovatti, MD |
| Sherrell Aston MD       | Jason Diamond, MD      | Ernest Layton, MD         | Aviva Preminger, MD    |
| Steven Bloch, MD        | Ronald Downs, MD       | Gilbert Lee, MD           | Mark Schusterman, MD   |
| Jason Bloom, MD         | Diane Duncan, MD       | Elie Levine, MD           | Anil Shah, MD          |
| Guillermo Blugerman, MD | Bryan Forley, MD       | Thomas Loeb, MD           | Spero Theodorou, MD    |
| Peter Bray, MD          | Yael Halaas, MD        | Constantino Mendieta, MD  | Glenn Vallecillos, MD  |
| Christopher Chia, MD    | Vu Ho, MD              | Brent Moelleken, MD       | Jennifer Walden, MD    |
| Michael Ciaravino, MD   | Alfredo Hoyos, MD      | James Newman, MD          | ...And more            |



**I N M O D E**  
*aesthetic solutions*

TOLL FREE {855} 411 2639 | [WWW.BODYTITE.COM](http://WWW.BODYTITE.COM)



DR. P. ROVATTI




DR. D. DUNCAN



DR. D. HURWITZ



DR. G. BLUGERMAN



BODYSCLUPT

DR. S. THEODOROU & DR. C. CHIA



## MEET THE BODYTITE PRO

### THE PREMIUM 40W EXPERIENCE

If you have patients that want exceptional contouring and lipolysis results without the excision, large scars, or downtime, then the BodyTitePRO is a great investment for your practice.

A key concern with the use of today's fat extraction technology is the compromised skin remaining in the treatment area. BodyTite's unprecedented three-dimensional contraction addresses those concerns.

Previously, patients only had the following options: multiple non-invasive treatment options with limited skin contraction results, a surgical tummy tuck or facelift with downtime, or a brachioplasty leaving large excision scars.

Patients who wanted less invasive procedures were left without a solution. With the BodyTite system, you can provide minimally invasive treatment solutions and bridge the gap to attract a whole new segment of patients to your practice.